

Molluscum Extractor

Description:

The **Molluscum Extractor** is a specialized surgical instrument used for the **removal of molluscum contagiosum lesions** and other skin-related dermatological extractions. Its unique round and oval loop design facilitates precise and safe extraction with minimal trauma to surrounding skin.

Specifications:

- **Sizes Available:**
 - Round/Oval loops in sizes 2/3, 3/4, and 4/5
- **Loop Shape:** Round and oval for various lesion sizes
- **Handle:** Ergonomic, textured for secure grip
- **Surface Finish:** Polished stainless steel for smooth contact

Sizes & Shapes:

- **Loop sizes:** 2/3, 3/4, 4/5 (mm approx. dimensions for lesion extraction)
- Available in **Round** and **Oval** loop shapes for flexibility during procedures

Material:

- Premium **surgical-grade stainless steel**
- Corrosion-resistant and autoclavable
- Durable and sterilizable for repeated use

Category:

- Surgical Dermatological Instruments

Product Form:

- Reusable instrument
- Non-disposable
- Available with loop variations as per sizes

Usage:

- Extraction of molluscum contagiosum lesions
- Dermatological procedures requiring lesion removal
- Suitable for both veterinary and human skin treatments

Advantages:

- Precise lesion removal with minimal skin trauma
- Ergonomic handle improves control and reduces hand fatigue
- Multiple loop sizes for different lesion dimensions

Disadvantages:

- Requires careful handling to avoid skin damage
- Must be properly sterilized before each use

Precautions:

- Sterilize before and after use
- Handle gently to prevent injury
- Avoid using excessive pressure during extraction

HS / HSN Code:

- **HS Code:** 9018.90
- **HSN Code:** 9018 (Surgical Instruments)

Handling Instructions:

- Hold firmly with a secure grip
- Use appropriate loop size depending on lesion
- Clean and sterilize thoroughly after use

Sterilization Details:

- Fully autoclavable/ Compatible with steam sterilization and chemical disinfectants

Veterinary Application:

Used in veterinary dermatology to remove skin lesions, growths, and other skin abnormalities in animals safely and effectively.

Human Application:

Utilized in dermatology clinics for removal of molluscum contagiosum and other skin lesions in patients.

Frequently Asked Questions (FAQs):

1. **What is a molluscum extractor used for?**
To safely remove molluscum contagiosum lesions and other skin abnormalities.
2. **What sizes are available?**
Round and oval loop sizes 2/3, 3/4, and 4/5.
3. **Can it be used on animals?**
Yes, it's suitable for veterinary dermatological procedures.
4. **How is it sterilized?**
It can be fully sterilized using steam autoclave or chemical disinfectants.
5. **Is the instrument reusable?**
Yes, made from surgical-grade stainless steel for repeated use.

