

Oral Feeding Handle Bottle – with Nipple (Silicone / Rubber)

Description

The **Oral Feeding Handle Bottle** is a veterinary-grade feeding device used to provide milk, nutrition, or medications to young or weak animals. It features an ergonomic handle for comfort and control. The bottle comes with either a **Silicone or Rubber nipple**, depending on the model selected, to cater to different animal sensitivities and feeding preferences.

Specifications

- **Material (Bottle):** Medical-grade plastic or stainless steel
- **Nipple Types Available:**
 - **Silicone Nipple** – Soft, transparent, durable, ideal for sensitive animals
 - **Rubber Nipple** – Traditional, slightly firmer, widely used
- **Bottle Sizes:**
 - 100 ml
 - 250 ml
 - 500 ml
 - 1000 ml
- **Handle:** Built-in side handle for grip and ease of use
- **Cap:** Secure screw-on cap for leak-proof feeding

Types

- Oral Feeding Bottle **with Silicone Nipple**
- Oral Feeding Bottle **with Rubber Nipple**
(Each type sold separately)

Product Form

- Bottle + Single Nipple (either Silicone or Rubber, based on selection)

Material

- **Bottle:** Polycarbonate or Stainless Steel
- **Nipples:** Medical-grade Silicone / Rubber

Category

- Veterinary Feeding Accessories

Veterinary Application

- Used for feeding milk, oral solutions, or supplements to neonates and sick animals (puppies, kittens, calves, lambs, etc.)

Human Application

- **Not Applicable**

Advantages

- **Two Nipple Options:** Choose based on animal's need and comfort
- **Easy to Use:** Ergonomic design for one-hand operation
- **Safe & Durable:** Made from veterinary-approved materials
- **Customizable Flow:** Nipple tips can be trimmed as needed

Disadvantages

- **Not Interchangeable in a Single Pack:** Must choose nipple type when ordering
- **Requires Regular Cleaning:** For hygiene and longevity

Precautions

- Always check and clean nipples before and after use
- Select the right nipple type for your animal's age and preference
- Replace worn-out nipples to avoid feeding issues

HS / HSN Code

- **HS Code:** 9018.39
- **HSN Code:** 9018.39.00
- **Description:** Veterinary surgical and feeding instruments

Sterilization Details

- **Silicone Nipples:** Boilable, disinfectant-safe
- **Rubber Nipples:** Clean with warm water, avoid over-boiling
- **Bottle:** Can be autoclaved (if stainless steel) or cleaned with hot water (if plastic)

Frequently Asked Questions (FAQs)

1. **Can I buy both silicone and rubber nipples?**
Yes, but they are sold as separate variants — not in one package.
2. **Which nipple is better?**
Silicone is softer and longer-lasting; rubber is firmer and more traditional.
3. **Are the bottles reusable?**
Yes, both the bottle and nipples are reusable with proper sterilization.
4. **What animals can this be used for?**
Commonly used for puppies, kittens, calves, piglets, and other neonates.
5. **Can this be used for humans?**
No, it is intended strictly for veterinary use.