

## Post Mortem Set – Bag

### Description

The **Post Mortem Set – Bag** is a comprehensive kit containing essential instruments required for conducting thorough post mortem examinations on small and large animals. Designed for veterinarians, veterinary pathologists, and diagnostic labs, these sets come in two variants — **18 No.** for small animals and **22 No.** for large animals — to cater to the size and complexity of different species.

### Specifications

- **Set Variants:**
  - **18 No. Set:** Instruments optimized for small animals (dogs, cats, rabbits, etc.)
  - **22 No. Set:** Instruments designed for large animals (cattle, horses, buffalo, etc.)
- **Material:** High-quality surgical stainless steel (SS 304 / 316)
- **Packaging:** Durable, portable bag for easy transport and storage
- **Contents:**
  - Scalpels and blades
  - Scissors (curved & straight)
  - Forceps (locking & non-locking)
  - Bone saw
  - Probes and hooks
  - Needles and suture materials
  - Other essential dissection tools
- **Bag Material:** Heavy-duty nylon or canvas with compartments for each instrument

### Product Form

- Complete Post Mortem Set in a portable carrying bag
- Individually sterilized instruments within the set
- Available in 18 No. and 22 No. configurations

## Veterinary Application

- Used for post mortem examinations of **small animals (18 No. set)** and **large animals (22 No. set)**
- Helps in diagnosis of cause of death, pathological studies, and research
- Suitable for veterinary clinics, diagnostic labs, and educational institutions

## Human Application

- **Not Applicable**

## Advantages

- Comprehensive sets covering all necessary post mortem instruments
- High-quality stainless steel ensures durability and corrosion resistance
- Portable and organized carrying bag for ease of use in the field
- Instruments designed for precision and safety during dissection

## Disadvantages

- Requires regular maintenance and sterilization
- Some instruments may require replacement over time due to wear

## Precautions

- Sterilize instruments before and after use
- Handle sharp instruments with care to avoid injury
- Store in a dry, clean environment to prevent rust
- Replace damaged instruments immediately

## HS / HSN Code

- **HS Code:** 9018.39
- **HSN Code:** 9018.39.10
- **Description:** Surgical instruments and appliances for veterinary use

## **Sterilization Details**

- Autoclave sterilization recommended for all instruments
- Clean thoroughly to remove tissue debris before sterilization
- Use protective gloves during handling and cleaning

## **Frequently Asked Questions (FAQs)**

1. **What's the difference between the 18 No. and 22 No. Post Mortem Sets?**  
The 18 No. set is tailored for small animals with fewer instruments, while the 22 No. set includes additional tools suitable for large animal dissection.
2. **Are the instruments reusable?**  
Yes, all instruments are reusable with proper cleaning and sterilization.
3. **What type of bag is used for the set?**  
Durable nylon or canvas bag with compartments for organized storage and easy transport.
4. **Can these sets be customized?**  
Customization is possible upon request to suit specific veterinary needs.
5. **Is training required to use all instruments in the set?**  
Basic veterinary knowledge is essential; specialized training is recommended for detailed post mortem procedures.

