

PPIUCD Tray Without Lid

Description

The **PPIUCD Tray without Lid** is a sturdy, sterilizable stainless steel tray designed to hold and organize instruments during postpartum intrauterine contraceptive device (PPIUCD) insertions. It offers a reliable sterile platform for clinical use, especially in gynaecological procedures.

Specifications

- **Material:** Medical-grade **Stainless Steel (SS 304 / SS 410)**
- **Tray Type:** Open stainless steel tray **without lid**
- **Surface Finish:** Matte / Mirror polish (as needed)
- **Standard Dimensions** (can vary):
 - 8 x 6 x 2.5 inches
 - 10 x 8 x 2.5 inches
- **Instrument Holders:** Moulded compartments or open layout
- **Reusability:** Fully autoclavable and reusable

Product Form

- Available as **empty tray only** or **custom pre-fitted set** (on request)
- No lid included

Veterinary Application

Not Applicable for veterinary use

Human Application

- Used in **postpartum IUCD insertions**
- Keeps instruments organized and sterile during procedures
- Suitable for hospitals, clinics, and labour rooms

✓ Advantages

- Durable stainless steel ensures longevity
- Easy to sterilize using autoclave or chemical methods
- Lightweight compared to trays with lids
- Customizable compartments available

✗ Disadvantages

- No lid means less protection from contamination during storage/transport
- Must be covered or handled carefully to maintain sterility

⚠ Precautions

- Always sterilize before and after use
- Store in a clean, dry place
- Handle with care to avoid contamination
- Use sterile gloves during handling

🔑 HS / HSN Code

- **HS Code:** 9018.90
- **HSN Code:** 9018.90.99

🧼 Sterilization Details

- Compatible with **autoclave (121–134°C)**
- Can be sterilized with ETO or chemical disinfectants
- Avoid abrasive scrubbing to maintain surface integrity

? Frequently Asked Questions (FAQs)

- 1. Is this tray supplied with instruments?**
No, usually sold empty or customized on request.
- 2. Can this tray be used without a lid safely?**
Yes, but it requires careful handling to maintain sterility.
- 3. What sizes does this tray come in?**
Common sizes: 8x6x2.5", 10x8x2.5" — custom sizes available.
- 4. Is this tray reusable?**
Yes, designed for repeated sterilization and reuse.
- 5. Is this tray applicable for veterinary procedures?**
No, it is intended for human gynaecological use only.

