

Milker Forceps

Description:

Milker Forceps are specialized veterinary instruments used during the **milking process or udder treatment** to handle the teats safely without causing injury or discomfort. They assist in managing teat placement, applying milking tubes, or controlling flow when treating mastitis or administering medication.

Specifications:

- Long-handled, smooth-tipped forceps
- Non-slip gripping surface
- Tension-adjustable mechanism (in some variants)
- Designed for precision without causing trauma to teat tissue

Sizes & Shapes:

- Standard size (approx. 6"–8")
- Straight or slightly curved tip
- Smooth or rubber-coated jaws (optional)

Material:

- **Stainless steel** body
- Optional rubber/silicone inserts for soft grip
- Rust and corrosion-resistant

Category:

- Veterinary Instruments
- Sub-category: Milking & Udder Care Tools

Product Form:

- Reusable
- Non-disposable
- Autoclavable

Usage:

- Handling and managing teats during milking or treatment
- Applying teat cannulas, plugs, or medication
- Supports mastitis treatment and teat care procedures
- Useful in dairy farms and vet practices

Advantages:

- Prevents direct hand contact with infected teats
- Reduces contamination risk
- Gentle grip to avoid damage or pain
- Durable and long-lasting tool

Disadvantages:

- Mishandling can still cause teat injury
- Must be kept sterile to avoid infections

Precautions:

- Always sterilize before use
- Avoid excessive pressure during application
- Use with gloves in clinical settings

HS / HSN Code:

- **HS Code:** 9018.90
- **HSN Code:** 9018 (Veterinary Instruments)

Handling Instructions:

- Hold forceps at the finger grips
- Use a steady hand for inserting or adjusting any teat-related tools
- Clean and dry thoroughly before storage

Sterilization Details:

- Fully autoclavable
- Can also be disinfected using veterinary-grade solutions

Veterinary Application:

Primarily used in **dairy veterinary care** to assist in mastitis treatment, milking aid, and gentle handling of sensitive udder and teat tissues.

Human Application:

 **Not applicable** – Designed exclusively for veterinary use.

Frequently Asked Questions (FAQs)

1. **What are Milker Forceps used for?**
For safe handling of animal teats during milking or treatment procedures.
2. **Are these forceps reusable?**
Yes, made from high-quality stainless steel and autoclavable.
3. **Can Milker Forceps prevent mastitis?**
Indirectly — by reducing contamination and handling trauma, they support udder health.
4. **Can this be used for goats and buffaloes too?**
Yes, suitable for all dairy animals including cows, goats, and buffaloes.
5. **Do they come in different sizes?**
Most come in standard sizes but variants may include curved or rubber-coated tips.