

Dressing Scissor

Description

Dressing Scissors are precision cutting instruments designed for veterinary and medical use. Available in both straight and curved blade options, these scissors are ideal for cutting bandages, dressings, and delicate tissues without causing damage. Their ergonomic design ensures control and accuracy during procedures.

Specifications

- **Blade Type:** Available in **Straight** and **Curved**
- **Sizes:** 4", 5", 6", 8", 10" (length in inches)
- **Material:** High-quality stainless steel, corrosion-resistant
- **Handle:** Ring handle design for secure grip
- **Finish:** Polished or matte finish
- **Tip:** Blunt tip or sharp tip options available

Usage

- Cutting bandages, dressings, and surgical tapes
- Trimming wounds and removing dead tissue gently
- Used in both veterinary and human medical applications
- Curved scissors offer better access in tight or awkward areas

Advantages

- **Versatile:** Suitable for various cutting tasks in medical and veterinary settings
- **Durable:** Made from stainless steel for longevity and resistance to corrosion
- **Precise:** Sharp blades ensure clean cuts without tearing tissues or dressings
- **Ergonomic:** Comfortable handles reduce hand fatigue during use

Disadvantages

- Requires regular sharpening to maintain blade efficiency
- Can cause injury if mishandled due to sharp edges
- Not suitable for cutting very tough materials

Precautions

- Always sterilize before and after use
- Handle with care to avoid accidental cuts
- Store in a safe place away from children
- Use only for intended purposes to prolong tool life

Handling Instructions

- Clean and sterilize with appropriate disinfectants after use
- Dry thoroughly to prevent rusting
- Sharpen blades regularly to maintain cutting quality
- Inspect before each use for damage or dullness

Sterilization Details

- Autoclaving recommended for sterilization
- Can be disinfected using chemical sterilant approved for surgical instruments
- Ensure thorough drying after sterilization to prevent corrosion

HSN & HS Code

- **HSN Code:** 9018.31
- **HS Code:** 9018.31.90
Surgical scissors for medical and veterinary use

Human Application:

Widely used in hospitals, clinics, and surgical centers for cutting dressings and tissues

Veterinary Application:

Used in veterinary clinics for dressing wounds, trimming tissues, and surgical procedures

? Frequently Asked Questions (FAQs)

1. What is the difference between straight and curved dressing scissors?

Straight scissors are ideal for general cutting, while curved scissors allow better access to hard-to-reach areas and provide precision in curved cuts.

2. Which size dressing scissors should I choose?

Size depends on the application—smaller sizes (4-6 inches) are for delicate work; larger sizes (8-10 inches) are for bigger animals or tougher dressings.

3. How do I sterilize dressing scissors?

They can be sterilized using an autoclave or chemical disinfectants suitable for surgical instruments.

4. Can dressing scissors be used for human and veterinary purposes?

Yes, these scissors are designed for both human medical and veterinary use.

5. How often should I sharpen the scissors?

Sharpen blades as soon as you notice any dullness affecting the cutting quality to maintain efficiency.

