

SS BOWL (CAPACITY TYPE)

DESCRIPTION


An SS Bowl (Stainless Steel Capacity Bowl) is a durable, corrosion-resistant medical utility bowl used for holding solutions, instruments, dressings, specimens, or cleaning fluids during surgical, medical, and veterinary procedures. Manufactured from high-grade stainless steel, it features a seamless design with a smooth, hygienic surface that prevents contamination and ensures easy cleaning. These bowls come in various capacities to suit different procedural needs in clinics, hospitals, and laboratories.

SPECIFICATION

- **Material:** High-grade stainless steel (SS 304 / SS 316)
- **Design:** Seamless round bowl with rolled rim
- **Capacity Range:** 50 ml to 2000 ml
- **Finish:** Polished or matte
- **Thickness:** 0.5 mm – 1.0 mm
- **Stability:** Flat base for secure placement
- **Usage Type:** Reusable, autoclavable
- **Hygiene Standard:** Rust-free, non-reactive surface

SIZES (CAPACITY OPTIONS)

- 50 ml
- 100 ml
- 150 ml
- 250 ml
- 350 ml
- 500 ml
- 750 ml
- 1 Litre
- 1.5 Litre
- 2 Litre



◆ SHAPES

- Standard round SS bowl
- Deep bowl
- Wide-mouth bowl
- Stackable nested bowls

🧪 MATERIAL TYPES

- Stainless Steel SS 304
- Stainless Steel SS 316 (premium surgical grade)
- High-polish surgical steel

📁 CATEGORY

- Medical Accessories
- Surgical & Instrument Bowls
- Veterinary Utility Equipment
- Hospital & Laboratory Containers

📄 PRODUCT FORM

- Single-piece stainless-steel bowl
- Non-sterile when supplied
- Reusable & easy to disinfect

🔧 USAGE

Human Application:

- Holding antiseptics, saline, cleaning solutions
- Collecting samples or small waste
- Storing instruments during procedures
- Dressing & wound-care setup
- Mixing solutions



Veterinary Application:

- Holding instruments & medications
- Mixing and preparing fluids
- Cleaning, soaking, and treatment support
- Cage-side and surgical table use

★ ADVANTAGES

- Rust-proof stainless steel
- Multiple capacity options
- Fully autoclavable
- Hygienic seamless design
- Durable & long-lasting
- Easy to wash and maintain
- Lightweight and stable

⚠ DISADVANTAGES


- Becomes hot after autoclaving
- May dent if dropped
- Not suitable for harsh corrosive chemicals
- Conducts heat quickly during procedures

🚫 PRECAUTIONS

- Avoid dropping to prevent dents
- Ensure bowl is properly dried to prevent water marks
- Use only for appropriate medical solutions
- Do not expose to extremely corrosive substances

📖 HANDLING DETAILS

- Hold by rolled rim for stability
- Cover when holding sterile solutions
- Stack carefully to prevent scratches



STERILISATION DETAILS

- Fully autoclavable
- EO (Ethylene Oxide) compatible
- Chemical disinfectant safe
- Ensure complete drying after sterilisation
- Avoid metal-scratching scrubbers

HS CODE (Harmonized System Code)

HS Code: 7323 – Stainless-steel household/medical articles.

HSN CODE (India)

HSN Code: 7323.93 – Stainless steel medical/hospital utility items.

FAQ (Frequently Asked Questions)

1. What is an SS Bowl used for?

It is used to hold solutions, dressings, instruments, or samples during medical and veterinary procedures.

2. Is it reusable?

Yes, it is fully reusable and autoclavable.

3. What capacities are available?

Multiple sizes from 50 ml to 2 litres.

4. Is it rust-proof?

Yes, made from high-grade stainless steel.

5. Can it be used for veterinary procedures?

Yes, widely used in all animal-care setups.