

# Pozzi Forceps

## Description:

**Pozzi Forceps**, also known as **traction forceps**, are specialized surgical instruments used to grasp and manipulate tough tissues such as fascia or uterine tissue. The 10-inch size offers extended reach and control during gynaecological and general surgical procedures.

## Specifications:

- Length: **10 inches**
- Curved design with three sharp, triangular prongs at the tip
- Robust handle for firm grip and precise control

## Sizes & Shapes:

- Available in **10-inch** length (standard for this type)
- Curved, serrated jaws with sharp prongs for secure tissue grip

## Material:

- Made from **high-quality surgical stainless steel**
- Corrosion-resistant and autoclavable

## Category:

- Surgical Instruments
- Sub-category: Tissue Grasping and Traction Instruments

## Product Form:

- Single-piece stainless steel forceps
- Reusable and sterilizable

### **Usage:**

- Used for grasping tough tissues during surgery, especially in gynaecology and obstetrics
- Provides strong traction on fascia, uterine tissue, or other dense tissues

### **Advantages:**

- Strong grip with minimal tissue damage
- Durable and long-lasting
- Easy to handle and manoeuvre during procedures

### **Disadvantages:**

- Sharp prongs can cause tissue trauma if not used carefully
- Requires skilled handling to avoid excessive damage

### **Precautions:**

- Use gently to avoid tearing tissues
- Ensure sterilization before use
- Regularly inspect for wear or damage

### **HS / HSN Code:**

- **HS Code:** 9018.39.00
- **HSN Code:** 9018 (Surgical instruments)

### **Handling Instructions:**

- Handle with sterile gloves
- Use with care to maintain tissue integrity
- Clean and autoclave after each use

### **Sterilization Details:**

- Compatible with autoclaving (134°C)
- Can also be sterilized using chemical disinfectants

### **Veterinary Application:**

Used in veterinary surgeries for tissue traction and manipulation, particularly in large animal gynaecological procedures.

### **Human Application:**

Widely used in human gynaecological and obstetric surgeries for tissue handling and traction.

### **Frequently Asked Questions (FAQs):**

1. **Is the Pozzi Forceps reusable?**  
Yes, made from surgical stainless steel and fully sterilizable.
2. **What tissues can it grasp?**  
Tough tissues like fascia, uterine tissue, and dense connective tissue.
3. **Can it be used in veterinary surgeries?**  
Yes, commonly used in large animal gynaecological procedures.
4. **Is the 10" size standard?**  
Yes, 10 inches is the typical length for Pozzi Forceps.
5. **How should it be sterilized?**  
Autoclaving at 134°C is recommended.